gneovo

MD-2702 CLINICAL REVIEW MONITOR




FEATURES

- Full HD 1920 x 1080 resolution
- 178° /178° wide viewing angle
- Support DICOM Part 14 with calibrated greyscale and brightness
- Crisp 20,000,000:1 dynamic contrast ratio .
- 3-sided narrow bezel design for seamless multi-display setups
- Medical adapter certified to the IEC/EN 60601-1 medical standard .
- White exterior design for easy contaminant detection
- Versatile Connectivity: DisplayPort, HDMI and VGA .
- . HDMI inputs for Full HD connectivity with a variety of devices
- Selectable aspect ratio .
- . VESA mounting for flexible installation in healthcare settings
- · Height adjustable stand allows users to tilt, pivot and swivel



Featuring DICOM preset and medical-grade power supply, the MD-2702 clinical review monitor ensures correct and consistent brightness and greyscale levels of medical images in detail at any viewing angle. Its image stability, medical-grade power supply, and VESAcompatible mounts give clinical professionals confidence in clinical review and easy to set up a mobile workstation while meeting the standards of premium quality assurance and financial benefits.

TECHNICAL DRAWING





SPECIFICATIONS

Panel		
Panel Type	LED-Backlit TFT LCD (VA Technology)	
Panel Size	27.0"	
Max. Resolution	FHD 1920 x 1080	
Pixel Pitch	0.311 mm	
Brightness	300 cd/m ²	
Lighthood	180 cd/m ² (DICOM Preset)	
Contrast Ratio	20,000,000:1 (DCR)	(Trose
Viewing Angle (H/V)	178°/178°	
Display Colour	16.7M) =) (
Response Time	5 ms	
Surface Treatment	Anti-Glare Treatment (Haze25%), 3H Hard Coating	
Frequency		
Frequency (H)	30 kHz-83 KHz	
Frequency (V)	50 Hz-76 Hz	
Signal Input		
DisplayPort	x 1	-
HDMI	1.4 x 1	
VGA	15-Pin D-Sub x 1	0 0
Power		
Power Supply	External	
Power Requirements	DC 12V, 2.5A	
On Mode	16W (On)	
Stand-by Mode	< 0.3W	
Off Mode	< 0.2W	
Operating Conditions		
Temperature	0°C-40°C (32°F-104°F)	
Humidity	10%-85% (non-condensing)	
Storage Conditions		
Temperature	'-25°C-55°C (-13°F-131°F)	
Humidity	5%-93%(non-condensing)	
Mounting		
VESA FPMPMI	Yes (100 x 100 mm)	
Stand		
Tilt	-5° to 23°	
Pivot	-2.5° to 92.5°	
Swivel	± 65°	
Height Adjustment	0-110 mm	
Security		
Kensington Security Slot	Yes	
Dimensions		
Product with Base (W x H x D)	612.6 x 519.8 x 249.0 mm (24.1" x 20.5" x 9.8")	
Product w/o Base (W x H x D)	612.6 x 367.1 x 50.4 mm (24.1" x 14.6" x 2.0")	
Packaging (W x H x D)	700.0 x 456.0 x 282.0 mm (27.6" x 18.0" x 11.1")	
Weight		
Product w/o Base	3.7 kg (8.2 lb)	
Product with Base	6.1 kg (14.5 lb)	
Packaging	8.8 kg (19.4 lb)	
Regulation Approval		
Certifications & Compliance	CE, LVD, WEEE, RoHS, REACH, EAC, UKCA, BSMI, FCC, ICES-003	
Power Adaptor Certifications	ADAPTER/ATM065T-P120: CE, FCC, PSE, WEEE, BSMI, Medical Grade	
Accessories		
Supplied	Power Adaptor, Power Cord, HDMI Cable, Quick Start Guide	
Optional Accessories	r ower Adaptor, i ower oord, ribiwi oabie, Quick Start duide	
	CMP-01+WMK-03	
Ceiling Mount (Landscape only)		
Wall Mount (Landscape only)	WMA-01, WMK-01, WMK-03, PMK-01	
Desk Mount (Landscape only)	DMC-01, ES-02	



All specifications are subject to change without prior notice. Copyright © 2023 AG Neovo. All rights reserved. The name AG Neovo is a trademark of AlC. All other trademarks are the property of their respective owners. P/C: MD272A Version 2.2, 20230220